

## High Speed Analysis of Estradiol in accordance with chapter 621 in USP 39

No. SCA\_190\_032

High throughput analysis has been advanced dramatically in recent years with the increasing necessity to improve productivity and operational efficiency. Especially HPLC has also been in the spotlight thanks to significant advances in ultra-high-speed analysis technology, in particular ultra-high performance LC and micro-particle column packing material. The recently revised General Chapter 621 of the United States Pharmacopoeia (USP 621) now permits a degree of adjustment of HPLC and GC parameters, specifically aimed to satisfy the requirements of system suitability. Taken account of USP 621, this Application News introduces an example of isocratic analysis of Estradiol monography in accordance with USP and still fulfilling the allowable adjustment criteria. Estradiol is a steroid and the primary female sex hormone. Additionally, an example of analysis that can be completed in a significantly shorter time than that described in the USP General Chapter 621 Chromatography is presented here.

### ■ Allowable adjustments to HPLC parameters

Table 1 shows the parameters which may be changed according to USP 621. The analysis was performed under isocratic conditions. Additionally, the actual permissible ranges within these LC parameters are shown.

### ■ Speed enhancement for USP method

The permissible ranges within the analytical conditions may be modified and are specified in the USP General Chapters: <621> Chromatography. Changing these analytical conditions within the range makes it possible to shorten the analysis time. For details regarding

changes that can be used to allow fast USP-compliant analysis, please refer to Application News L464. Shortening analysis time can be accomplished in three ways, 1) by shortening the column, 2) by lowering the inner diameter and 3) by increasing the flowrate while maintaining the linear velocity. To preserve the resolution of the separation, the column length and particle size may be modified as long as the ratio of L (column length) to dp (column particle size) remains in the specified range (permissible range: -25 % to +50 %). For the original USP method, a column with the dimensions 300 mmL. x 3.9 mm I.D., and 10 µm particle size was used. We selected a column size of 100 mmL. x 3.0 mm I.D., and 3 µm particle size with constant L/dp ratio (Table 2). For further details, please see Table 3. The flowrate, proportional to the column cross-sectional area, and inversely proportional to the particle diameter (see text for permissible limits), was determined as 1.5 mL/min. Table 4 shows the analytical conditions.

Table 1: Allowable adjustments to HPLC parameters according to USP 621

|                   |   |
|-------------------|---|
| Particle size(dp) | L/dp ratio constant or Theoretical plate number: -25 to + 50%     |
| Column length(L)  |   |
| Column ID(dc)     | Any allowed if linear velocity is constant                        |
| Flowrate          | Combination* of dp and dc : ±50%                                  |
| Injection Vol.    | Can be adjusted as consistent with precision and detection limits |
| Column Temp.      | ±10 °C  |

$$* F_2 = F_1 \times [(dc_2^2 \times dp_1)/(dc_1^2 \times dp_2)]$$

F<sub>1</sub> and F<sub>2</sub> are the flow rates for the original and modified conditions, respectively; dc<sub>1</sub> and dc<sub>2</sub> are the respective column diameters; and dp<sub>1</sub> and dp<sub>2</sub> are the particle sizes.

Table 2: Selection of column for speed enhancement

|                        | Column size           | L/dp  | Ratio          |
|------------------------|-----------------------|-------|----------------|
| USP<br>Original Method | 300 x 3.9 mm<br>10 µm | 30000 | 1<br>(100%)    |
| USP<br>Fast Method     | 100 x 3.0 mm<br>3 µm  | 33333 | 1.11<br>(+11%) |

Table 3: Column selection for speed enhancement in case of fixing the particle size and column I.D.

|                   | USP method of | Allowable range    | Modified method |
|-------------------|---------------|--------------------|-----------------|
| Particle size(dp) | 10 µm         | 3 µm               | 3 µm            |
| Column ID(dc)     | 3.9 mm        | 3 mm               | 3 mm            |
| Column length(L)  | 300 mm        | 67.5 - 135 mm      | 100 mm          |
| Flowrate          | 1.0 mL/min    | 0.99 – 2.96 mL/min | 1.5 mL/min      |
| Injection Vol.    | 20 µL         | Variable           | 10 µL           |
| Column Temp.      | Unspecified   | Variable           | 40 °C           |

Table 4: Analytical conditions

|                |  |
|----------------|--|
| System         | LC-2040C 3D  |
| Column         | (1) Shim-pack GIS C18 (100 x 3.0 mm, 3 µm) (2) Shim-pack GIS C18 (250 x 4.0 mm, 10 µm) |
| Mobile phase   | A) Acetonitrile<br>B) Water A/B = 55/45 (v/v)  |
| Flowrate       | (1) 1.5 mL/min; (2) 1.0 mL/min   |
| Column Temp.   | 40 °C  |
| Injection Vol. | (1) 10 µL; (2) 20 µL   |
| Detection      | LC-2030/2040 PDA; D2 at 190-350  |

## ■ Results

The results are shown in Figure 1 and 2 and in Table 5. The speed enhancement is shown in Figure 1. Here, the retention times are much shorter with 0.7 minutes for Ethylparabene, 0.9 minutes for Estradiol and 1.2 minutes for Estrone, compared to the retention times in Figure 2 (4.1 minutes for Ethylparaben, 5.3 minutes for Estradiol and 7.0 minutes for Estrone).




Figure 1: Chromatogram of USP fast method for Estradiol (Peak 2) (20.01 µg/mL) with column (1) Shim-pack GIS C18 (100 x 3.0 mm; 3 µm). Peak 1: Ethylparaben (15.00 µg/mL); Peak 3: Estrone (12.11 µg/mL)




Figure 2: Chromatogram of USP original method for Estradiol (Peak 2) (20.01 µg/mL) with column (2) Shim-pack GIS C18 (250 x 4.0 mm; 10 µm). Peak 1: Ethylparaben (15.00 µg/mL); Peak 3: Estrone (12.11 µg/mL)

Table 5: Results of system suitability test using USP Method (original method and fast method).

| System Suitability       |                   | Ref. Value | USP original Value | Fast method Value |     |      |
|--------------------------|-------------------|------------|--------------------|-------------------|-----|------|
| Relative retention times | Estradiol         | 1.0        | 1.0                | 1.0               |     |      |
|                          | Estrone           | 0.7        | 0.8                | 0.8               |     |      |
|                          | Ethylparaben      | 1.3        | 1.3                | 1.3               |     |      |
| Resolution               | Estradiol-Estrone | ≥2.0       | 5.6                | 4.6               |     |      |
| RSD%                     | Estradiol         | ≤2%        | 6.2                | Pass              | 5.6 | Pass |

## ■ Conclusion

With the fast USP method, the original USP method, according to the reference value, was improved because the analysis time is shorter and solvent consumption is reduced. Ongoing with this, the cost per analysis is reduced significantly. Additionally, both column conditions are better than the USP reference.