

# Application Data Sheet

## No. 137

### GC-MS

Gas Chromatograph Mass Spectrometer

## Analysis of VOC and SVOC Emissions from Automotive Interior Materials in Accordance with VDA278 Using the Thermal Desorption Method

In recent years, measures to reduce the use of organic compounds in automotive interiors have progressed. In Germany, the VDA278 standards were created for the analysis of volatile organic compounds (VOC) and semi-volatile organic compounds (SVOC) produced from automotive interior materials. In the VDA278, measurement samples are added to a TD glass tube. The VOC (up to C20) and SVOC (up to C32) are heated at different temperatures, and the gases produced are loaded into a GC-MS. The VOC and SVOC from the automotive interior materials can be analyzed conveniently and quickly. However, since the gas is loaded directly, if the SVOC content is highly concentrated, caution is needed regarding carryover.

The TD-30 thermal desorption system features an inert sample line that is kept as short as possible, and samples are heated up to 300 °C, reducing SVOC carryover. In this investigation, an analysis of VOC and SVOC emissions from automotive interior materials was attempted in accordance with VDA278 using the TD-30.


### Experiment

Automotive interior materials (rubber, plastic, and leather) were sliced thinly, and TD glass tubes (from Shimadzu) were filled with approximately 30 mg of these samples. Both ends were fastened with 5 mg of quartz wool. The VOC samples were heated at 60 °C for 30 minutes, the SVOC samples were heated at 90 °C for 60 minutes, and the gases produced were loaded into a GC-MS. The samples loaded were analyzed in GC-MS Scan mode. The analysis conditions are shown in Table 1, and the analysis samples are shown in Fig. 1.

Table 1: Analysis Conditions

[Instrument Configuration]			
GC-MS:	GCMS-QP™ 2020		
Sample Loader:	TD-30R		
Workstation (GCMS-QP™2020):	GCMSsolution™ Ver.4.45		
Workstation (TD-30):	TD-30 Control Software		
Column:	SH-Rxi™-5Sil MS (60 m x 0.25 mm I.D., df = 0.25 μm) (SHIMADZU)		
[TD-30]		[GC]	
Tube Desorption Temperature:	90 °C for 30 min (VOC) 120 °C for 60 min (SVOC)	Control Mode:	Pressure
Tube Desorption Flowrate:	60 mL/min	Pressure:	200 kPa
Trap Cooling Temperature:	-20 °C	Injection Mode:	Split 1:100 (Column Flowrate 1.99 mL/min)
Trap Desorption Temperature:	280 °C for 10 min	Column Oven Temperature:	40 °C (3 min) – (10 °C/min) – 300 °C (13.5 min)
Joint Temperature:	280 °C	[MS]	
Valve Temperature:	250 °C	Ion Source Temperature:	200 °C
Transfer Line Temperature:	280 °C	Interface Temperature:	250 °C
		Measurement Mode:	Scan Measurement
		Scan Mass Range:	m/z 35-400
		Scan Event Time:	0.5 sec
		Scan Speed:	769 u/sec

#### (1) Flat Sample Tray Capable of Heating the Sample Directly



#### (2) Inert Sample Line as Short as Possible and Capable of Being Heated



The TD-30 sample line is designed to be as short as possible, and is not connected to the GC sample vaporization chamber or other unnecessary parts. In addition, all parts of the sample line can be heated to at least 300 °C, so carryover is not a concern, even for SVOC with high boiling points.

Fig. 1: Features of Analysis Samples and the TD-30 When Analyzing Automotive Interior Materials

## Analysis Results

### Evaluative Results for Calibration Curves and Recovery Rates

The standard samples for the calibration curves were prepared by diluting toluene and n-hexadecane with methanol to concentrations of 0.5 µg/µL. 4 µL of the sample was added to a Tenax® TA collection tube, and the response factor was calculated. The response factor was used in the calculation of the quantitative values of the compounds in the automotive interior materials. The formula is shown below.

In addition, in order to evaluate the recovery rate for the analysis system, a standard sample of typical VOCs (with a concentration of approximately 0.11 µg/µL) was prepared. 4 µL of this was added to a Tenax® TA collection tube and then analyzed. When the recovery rate was calculated from the response factor, values between 60 % and 140 % were obtained regardless of the compound, which is a favorable recovery rate.

$$R_f = \frac{\mu\text{g Toluene (C16)} \times 1000000}{\text{Peak area}}$$

Formula 1: Formula for the Response Factor (Rf)

$$\text{Emission } [\mu\text{g/g}] = R_f (\text{Toluene, C16}) \times \frac{\text{Peak area [count]}}{1000 \times \text{sample weight [mg]}}$$

Formula 2: Formula for the Quantitative Values (Emission[µg/g]) of Compounds Produced by Automotive Interior Materials

Table 2: Recovery Rate for Typical VOCs

Name of Compound	Recovery Ratio (%)
Benzene	106.53
Toluene	93.49
p-Xylene	99.91
o-Xylene	75.38
2-Ethyl-1-hexanol	101.26
2,6-Dimethylphenol	94.19
Dicyclohexylamine	89.15

### Analysis Results for Automotive Interior Materials

The quantitative values (µg/g) for compounds produced by rubber, plastic, and leather samples are shown in Table 3. A high concentration of Bis(2-ethylhexyl) phthalate (at a concentration of 333.28 µg/g) was detected from the leather sample. When a blank sample was analyzed immediately after analyzing the leather sample, the carryover was less than 0.05 %, which is a favorable result.

Table 3: List of Quantitative Values of Compounds Produced by Automotive Interior Materials

Name of Compound	VOC			SVOC		
	Rubber	Plastic	Leather	Rubber	Plastic	Leather
C8	0.00	0.00	0.00	0.00	0.00	0.11
Toluene	0.35	0.54	0.53	0.31	0.44	0.24
C9	0.00	0.00	0.00	0.00	0.00	0.13
C11	0.00	0.00	0.00	0.00	0.00	0.31
Benzene, 1,3-dichloro-	0.00	0.00	0.00	0.00	0.00	0.08
2-Propyl-1-pentanol	0.36	0.52	0.73	0.11	0.18	0.78
C12	0.00	0.00	0.17	0.00	0.03	0.06
Nonanal	0.00	0.00	0.43	0.09	0.06	0.87
C13	0.20	0.14	0.26	0.09	0.13	0.13
C15	0.14	0.12	0.36	0.13	0.16	0.14
C16	0.31	0.00	0.60	0.42	0.16	0.86
C18	0.14	0.00	0.73	0.39	0.00	2.02
C19	0.00	0.00	0.30	0.39	0.00	1.37
Dibutyl phthalate	0.00	0.00	2.92	0.00	0.00	17.53
C20	0.00	0.00	0.18	0.14	0.00	1.28
C22	0.00	1.09	0.17	0.00	0.00	0.82
C23	0.00	0.00	0.15	0.00	0.00	0.82
C25	0.00	0.00	0.00	0.00	0.00	1.78
Bis(2-ethylhexyl) phthalate	0.41	1.60	33.67	0.00	0.00	333.28





Fig. 2: Chromatograms for Bis(2-ethylhexyl) Phthalate in the Leather and Blank Samples

GCMS-QP and GCMSsolution are trademarks of Shimadzu Corporation.  
Rxi is a registered trademark of Restek Corporation.  
Tenax is a registered trademark of Buchem B.V.

First Edition: December, 2017



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.